
DE 838 Diffractive Optical Element



Element number: DE 838

Product revision: A

Description: matrix 3 x 5 dots

Substrate material: fused silica

AR coating on both sides of the substrate: R < 0.5% within recommended wavelength range

Substrate size: 15.0 mm x 14.1 mm

Thickness: 2.3 mm

Design wavelength: 532 nm

Recommended wavelength range: 488 nm - 532 nm *

Typ. diffraction efficiency: 80% at design wavelength

Within the recommended wavelength range, the central spot / zeroth order (Z0) has a similar power as the desired offaxis orders of the dot matrix. Pattern size and pattern angles, and the ratio between central spot / zeroth order and desired orders will vary most with the wavelength. Diffraction efficiencies given on this datasheet have been measured using elements of product revision A.


The DOEs are best used with collimated or convergent laser sources. The microstructure surface should be oriented towards the laser. The 0-order spot is equivalent in size and shape to the original beam, but its power is attenuated.


Diffraction angles & efficiencies

Wavelength	Pattern Size @ 100 mm Distance				Pattern Angles			
λ [nm]	a [mm]	b [mm]	c [mm]	d [mm]	α [°]	β [°]	γ [°]	δ [°]
488	5.7	2.6	1.28	5.1	3.3	1.5	0.74	2.9
515	6.1	2.7	1.36	5.4	3.5	1.6	0.78	3.1
532	6.3	2.8	1.40	5.6	3.6	1.6	0.80	3.2

Table 1: Pattern size and pattern angle depending on the wavelength

Orders at 532nm






*the recommended wavelength range is defined with $4.3\% \le Z0 \le 6.3\%$ (mean off axis order ~5.3% ± 1%)



HOLOEYE Photonics AG Volmerstr. 1 12489 Berlin, Germany doe@holoeye.com www.holoeye.com



For testing or setups under

laboratory conditions, we offer

a version mounted in a black

anodized 25 mm aluminum frame for use with standard

Ø 25,00

25 mm anodized aluminum mount

with 14.0 x 13.1 mm clear aperture

The laser can be collimated for

long-range use or converging

for a fixed working distance.

size/thickness of each spot or

line depends on the focusing of

Please note that the

the laser.

laboratory holders.